

## Blunt-ended Ligation

*The maximum number of "correct" clones can generally be obtained from ligation reactions containing equimolar amounts of plasmid and target DNAs, with the total DNA concentration being <100 µg/ml. Blunt-end ligation catalyzed by bacteriophage T4 DNA ligase is suppressed by high concentrations (5 mM) of ATP and polyamines such as spermidine.*

## Buffers and Solutions

ATP (10 mM)

*Omit ATP from the ligation reaction in Step 5 if the ligation buffer contains ATP.*

PEG 8000 (30% w/v)

Sodium acetate (3 M, pH 5.2)

TE (pH 8.0)

1. In separate reactions, digest 1-10 µg of the plasmid DNA and insert DNA with the appropriate restriction enzyme(s) that generate blunt ends.
2. Purify the digested insert DNA and vector DNA by extraction with phenol:chloroform and ethanol precipitation.
3. Dissolve the precipitated DNAs separately in TE (pH 8.0) at a concentration of approx. 100 ng/µl. Calculate the concentration of the DNAs assuming that 1 bp has a mass of 650 daltons.
4. Dephosphorylate the plasmid vector DNA
5. Transfer appropriate amounts of the DNAs to microfuge tubes as follows:

vector\* (60 fmoles [approx. 100 ng])

insert (60 fmoles [approx. 10 ng])

linearized vector (contains 5'-terminal phosphates) (60 fmoles)

superhelical vector (6 fmoles [approx. 10 ng])

*\*Vector DNA is dephosphorylated*

*The molar ratio of plasmid vector to insert DNA fragment should be approx. 1:1 in the ligation reaction. The total DNA concentration in the ligation reaction should be approx. 10 ng/µl.*

add:

10x Ligation buffer 1.0 µl

Bacteriophage T4 DNA ligase 0.5 units

5 mM ATP 1.0 µl

H<sub>2</sub>O to 8.5 µl

30% PEG 8000 1-1.5 µl

*To achieve the maximum efficiency of ligation, set up the reactions in as small a volume as possible (5-10 µl).*

*The DNA fragments can be added to the tubes together with the H<sub>2</sub>O and then warmed to 45°C for 5 minutes to help dissociate any clumps of DNA that have formed during fragment preparation. Chill the DNA solution to 0°C before the remainder of the ligation reagents are added.*

6. Incubate the reaction mixtures overnight at 16°C or for 4 hours at 20°C.
7. Transform competent *E. coli* with the ligation reactions. As controls, include known amounts of a standard preparation of superhelical plasmid DNA to check the efficiency of transformation, vector without insert and vector without dephosphorylation.