

## Galactose induction for analytical purposes

- solutions:
  - 5x buffer A:
    - 0.1 M Tris pH 7.5
    - 5 mM EDTA
    - 50% glycerol
  - extraction buffer:
    - 1x buffer A
    - 150 mM NaCl
    - (0.2 mM PMSF [toxic!])
  - basic medium:
    - 0.67% yeast nitrogen base w/o amino acids 6.7 g
    - amino acid drop out mix 0.87 g
    - 2% (w/v) sodium lactate 25.6 ml 60%(w/w)
    - 3% glycerol 30 ml
    - add H<sub>2</sub>O to 1 l

- preculture in 5 ml SD medium, incubate 24 h at 30°C
- start main culture in 50 ml basic medium at OD<sub>600</sub>≈~0.2
- incubate ~ 16h at 30°C (OD<sub>600</sub> should reach 1-2)
- add 2% galactose (5 ml 20%)
- induce for 6 h at 30°C
- spin down at ~4000 rpm for 10 min
- resuspend pellet in 1 ml extraction buffer and transfer to screw cap tube for cell disruption
- spin down in table top centrifuge
- pellet can be frozen in liquid N<sub>2</sub> and stored at -80°C

## Protein extraction

- material:
  - acid washed glass beads (Ø=0.45 mm):
  - leave 15 min in 5 M nitric acid (HNO<sub>3</sub>)
  - wash with H<sub>2</sub>O until pH is neutral (check with pH paper)
- add 250 µl extraction buffer to cells (make sure it contains PMSF)
- add 250 µl acid washed glass beads
- extract proteins in bead beater (fast prep) for 1-2x 45 sec (setting 4)
- spin down 5-10 min at 4°C (ependorf centrifuge)
- transfer supernatant to fresh tube (~200 µl)
- measure protein concentration with Bradford assay (see Bio-Rad protein kit)
- store extracts at -80°C