

DNA strand exchange reaction

- solutions:
- 5x strand exchange buffer:
for 1 ml 5x SEB mix

150 μ l	1 M Tris Ac pH 7.5 [150 mM]
12.5 μ l	0.4 M DTT [5 mM]
25 μ l	10 μ g/ μ l BSA (from NEB) [250 μ g/ml]
125 μ l	0.1 M ATP [12.5 mM]
100 μ l	1 M MgOAc [100 mM]
100 μ l	1 M phosphocreatine [100 mM]
487.5 μ l	H ₂ O

 - stop buffer:
60 μ l 10% SDS
600 μ l 0.5 M EDTA
180 μ l 20 mg/ml Proteinase K
 - 10x loading buffer: bromphenolblue in 50% glycerol (no xylencyanol)
 - DNA
X174 ssDNA 130 ng/ μ l (412 μ M [ntd]) in TE
X174 dsDNA digested with *Pst*I 260 ng/ μ l
(412 μ M [bp]) in TE

Preparation of linear X174 dsDNA:

- digest 26 μ g of RFI dsDNA in 100 μ l of NEB buffer #3 with *Pst*I
- remove restriction enzyme with phenol/isoamyl alcohol extraction (2x),
precipitate with ethanol/NH₄Ac, dissolve in 50 μ l TE
- measure A₂₆₀ and determine concentration using absorption coefficient $\epsilon_{260}=6500/\text{cm}$
- adjust concentration to 412 μ M [bp] with TE

Reactions:

- mix:
2.2 μ l 0.1 M ATP
2.5 μ l 5x SEB
0.296 μ l 0.1 M spermidine
1.2 μ l 1 μ g/ μ l creatine kinase
1 μ l ssDNA
x μ l Rad51 [10.3 μ M]
add H₂O to 10.5 μ l
- incubate 15 min at 30°C
- add RPA [1.8 μ M]
- incubate 30 min at 30°C
- add 0.5 μ l dsDNA, Rad54/storage buffer [0.2 μ M] and 0.296 μ l 0.1 M spermidine
- incubate at 30°C for 1-6 h

- stop reaction by adding 2 μ l stop buffer
- incubate 20 min at 30°C
- add 2 μ l 10x loading buffer
- separate samples on 0.8% agarose gel (without ethidium bromide) at 25-30 V over night (100 ml gel)
- stain gel in 1 μ g/ml ethidium bromide for 1 h, destain in H₂O for 1-2 h,
- take photograph on UV transilluminator